

Product Datasheet

Human aMACR ELISA Kit (orb1199529)

Description	Human Alpha-Methylacyl Coenzyme A Racemase (aMACR) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UHK6
Sensitivity	0.059ng/mL
Entrez	23600
Expiration Date	Please enquire.