

Product Datasheet

Human ALDH1A1 ELISA Kit (orb1199538)

Description	Human Aldehyde Dehydrogenase 1 Family, Member A1 (ALDH1A1) ELISA Kit
Reactivity	Human
Range	0.625-40U/L
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P00352
Sensitivity	0.236U/L
Entrez	216
Expiration Date	Please enquire.