

Product Datasheet

Human AGR2 ELISA Kit (orb1199548)

Description	Human Anterior Gradient Protein 2 (AGR2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O95994
Sensitivity	0.061ng/mL
Entrez	10551
Expiration Date	Please enquire.