

## Product Datasheet

### Human ADRbK1 ELISA Kit (orb1199554)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Human Adrenergic Receptor Beta Kinase 1 (ADRbK1) ELISA Kit    |
| <b>Reactivity</b>          | Human                                                         |
| <b>Range</b>               | 0.312-20ng/mL                                                 |
| <b>Tested Applications</b> | ELISA                                                         |
| <b>Note</b>                | For research use only                                         |
| <b>Application notes</b>   | sandwich                                                      |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>P25098</b>                                                 |
| <b>Sensitivity</b>         | 0.105ng/mL                                                    |
| <b>Entrez</b>              | <b>156</b>                                                    |
| <b>Expiration Date</b>     | Please enquire.                                               |