

Product Datasheet

Human ADAMTS9 ELISA Kit (orb1199567)

Description	Human A Disintegrin And Metalloproteinase With Thrombospondin 9 (ADAMTS9) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9P2N4
Sensitivity	0.109ng/mL
Entrez	56999
Expiration Date	Please enquire.