

## Product Datasheet

### Human ADAM22 ELISA Kit (orb1199571)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Human A Disintegrin And Metalloprotease 22 (ADAM22) ELISA Kit |
| <b>Reactivity</b>          | Human                                                         |
| <b>Range</b>               | 0.156-10ng/mL                                                 |
| <b>Tested Applications</b> | ELISA                                                         |
| <b>Note</b>                | For research use only                                         |
| <b>Application notes</b>   | sandwich                                                      |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>Q9P0K1</b>                                                 |
| <b>Sensitivity</b>         | 0.057ng/mL                                                    |
| <b>Entrez</b>              | <b>53616</b>                                                  |
| <b>Expiration Date</b>     | Please enquire.                                               |