

Product Datasheet

Human ACAT2 ELISA Kit (orb1199580)

Description	Human Acetyl Coenzyme A Acetyltransferase 2 (ACAT2) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9BWD1
Sensitivity	31.5pg/mL
Entrez	39
Expiration Date	Please enquire.