

Product Datasheet

Human ABCC12 ELISA Kit (orb1199591)

Description	Human ATP Binding Cassette Transporter C12 (ABCC12) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q96J65
Sensitivity	0.072ng/mL
Entrez	94160
Expiration Date	Please enquire.