

Product Datasheet

Human AASDHPPT ELISA Kit (orb1199595)

Description	Human Aminoadipate Semialdehyde Phosphopantetheinyl Transferase (AASDHPPT) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9NRN7
Sensitivity	0.054ng/mL
Entrez	60496
Expiration Date	Please enquire.