

Product Datasheet

Porcine VEGFR1 ELISA Kit (orb1199606)

Description	Porcine Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	F1RSU5
Sensitivity	0.115ng/mL
Entrez	396763
Expiration Date	Please enquire.