

Product Datasheet

Rat SST ELISA Kit (orb1199696)

Description	Rat Somatostatin (SST) ELISA Kit
Reactivity	Rat
Range	7.81-500pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P60042
Sensitivity	2.16pg/mL
Entrez	24797
Expiration Date	Please enquire.