

Product Datasheet

Rat SPA ELISA Kit (orb1199700)

Description	Rat Surfactant Associated Protein A (SPA) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates or other biological fluids.
Uniprot ID	P08427
Sensitivity	6.9pg/mL
Entrez	24773
Expiration Date	Please enquire.