

Product Datasheet

Mouse SOD1 ELISA Kit (orb1199705)

Description	Mouse Superoxide Dismutase 1, Soluble (SOD1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P08228
Sensitivity	5.02pg/mL
Entrez	20655
Expiration Date	Please enquire.