

Product Datasheet

Rat SELL ELISA Kit (orb1199716)

Description	Rat Selectin, Leukocyte (SELL) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P30836
Sensitivity	0.041ng/mL
Entrez	29259
Expiration Date	Please enquire.