

Product Datasheet

Rat S100B ELISA Kit (orb1199725)

Description	Rat S100 Calcium Binding Protein B (S100B) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P04631
Sensitivity	12.1pg/mL
Entrez	25742
Expiration Date	Please enquire.