

Product Datasheet

Mouse PYY ELISA Kit (orb1199742)

Description	Mouse Peptide YY (PYY) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9EPS2
Sensitivity	4.9pg/mL
Entrez	217212
Expiration Date	Please enquire.