

Product Datasheet

Rat PROCR ELISA Kit (orb1199754)

Description	Rat Protein C Receptor, Endothelial (PROCR) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q4V8I1
Sensitivity	0.64ng/mL
Entrez	362248
Expiration Date	Please enquire.