

Product Datasheet

Rat PPARg ELISA Kit (orb1199763)

Description	Rat Peroxisome Proliferator Activated Receptor Gamma (PPARg) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O88275
Sensitivity	0.058ng/mL
Entrez	25664
Expiration Date	Please enquire.