

Product Datasheet

Rat NTXI ELISA Kit (orb1199817)

Description	Rat Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit
Reactivity	Rat
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	0.88ng/mL
Expiration Date	Please enquire.