


Product Datasheet

Rat NT-ProBNP ELISA Kit (orb1199819)

Description	Rat N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	27pg/mL
Expiration Date	Please enquire.



Standard
Curve of Rat
NT-ProBNP
ELISA Ki...