

Product Datasheet

Mouse MyD88 ELISA Kit (orb1199850)

Description	Mouse Myeloid Differentiation Factor 88 (MyD88) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P22366
Sensitivity	29pg/mL
Entrez	17874
Expiration Date	Please enquire.