

Product Datasheet

Rat MUC5AC ELISA Kit (orb1199852)

Description	Rat Mucin 5 Subtype AC (MUC5AC) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9ESP3
Sensitivity	0.115ng/mL
Entrez	682837
Expiration Date	Please enquire.