

Product Datasheet

Rat MMP1 ELISA Kit (orb1199875)

Description	Rat Matrix Metalloproteinase 1 (MMP1) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P81563
Sensitivity	6.7pg/mL
Entrez	300339
Expiration Date	Please enquire.