

Product Datasheet

Rabbit MIP1b ELISA Kit (orb1199886)

Description	Rabbit Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P46632
Sensitivity	6.2pg/mL
Entrez	100008984
Expiration Date	Please enquire.