

Product Datasheet

Rabbit INS ELISA Kit (orb1199948)

Description	Rabbit Insulin (INS) ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01311
Sensitivity	0.043ng/mL
Entrez	100009181
Expiration Date	Please enquire.