

Product Datasheet

Rat INHA ELISA Kit (orb1199955)

Description	Rat Inhibin A (INHA) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17490
Sensitivity	6.5pg/mL
Entrez	24504
Expiration Date	Please enquire.