

Product Datasheet

Rat GSTa3 ELISA Kit (orb1199990)

Description	Rat Glutathione S Transferase Alpha 3 (GSTa3) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P46418
Sensitivity	0.079ng/mL
Entrez	494500
Expiration Date	Please enquire.