

Product Datasheet

Mouse gGT1 ELISA Kit (orb1200011)

Description	Mouse Gamma-Glutamyltransferase 1 (gGT1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q60928
Sensitivity	0.056ng/mL
Entrez	14598
Expiration Date	Please enquire.