

Product Datasheet

Mouse FABP3 ELISA Kit (orb1200043)

Description	Mouse Fatty Acid Binding Protein 3, Muscle And Heart (FABP3) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11404
Sensitivity	26.3pg/mL
Entrez	14077
Expiration Date	Please enquire.