

## Product Datasheet

### Mouse CD4 ELISA Kit (orb1200131)

|                            |                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Cluster Of Differentiation 4 (CD4)<br>ELISA Kit                                 |
| <b>Reactivity</b>          | Mouse                                                                                 |
| <b>Range</b>               | 0.312-20ng/mL                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                 |
| <b>Note</b>                | For research use only                                                                 |
| <b>Application notes</b>   | sandwich                                                                              |
| <b>Sample Types</b>        | tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>P06332</b>                                                                         |
| <b>Sensitivity</b>         | 0.116ng/mL                                                                            |
| <b>Entrez</b>              | <b>12504</b>                                                                          |
| <b>Expiration Date</b>     | Please enquire.                                                                       |