

Product Datasheet

Mouse C5aR1 ELISA Kit (orb1200165)

Description	Mouse Complement Component 5a Receptor 1 (C5aR1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P30993
Sensitivity	0.128ng/mL
Entrez	12273
Expiration Date	Please enquire.