

Product Datasheet

Mouse ATM ELISA Kit (orb1200191)

Description	Mouse Ataxia Telangiectasia Mutated (ATM) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q62388
Sensitivity	5.6pg/mL
Entrez	11920
Expiration Date	Please enquire.