

Product Datasheet

Mouse APOA5 ELISA Kit (orb1200210)

Description	Mouse Apolipoprotein A5 (APOA5) ELISA Kit
Reactivity	Mouse
Range	3.12-200pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8C7G5
Sensitivity	1.45pg/mL
Entrez	66113
Expiration Date	Please enquire.