

Product Datasheet

Mouse APOA4 ELISA Kit (orb1200211)

Description	Mouse Apolipoprotein A4 (APOA4) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06728
Sensitivity	38pg/mL
Entrez	11808
Expiration Date	Please enquire.