

## Product Datasheet

### Mouse ANGPTL4 ELISA Kit (orb1200219)

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Angiopoietin Like Protein 4 (ANGPTL4) ELISA Kit                                                |
| <b>Reactivity</b>          | Mouse                                                                                                |
| <b>Range</b>               | 0.156-10ng/mL                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                |
| <b>Note</b>                | For research use only                                                                                |
| <b>Application notes</b>   | sandwich                                                                                             |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>Q9Z1P8</b>                                                                                        |
| <b>Sensitivity</b>         | 0.065ng/mL                                                                                           |
| <b>Entrez</b>              | <b>57875</b>                                                                                         |
| <b>Expiration Date</b>     | Please enquire.                                                                                      |