

Product Datasheet

Mouse ANGPTL2 ELISA Kit (orb1200220)

Description	Mouse Angiopoietin Like Protein 2 (ANGPTL2) ELISA Kit
Reactivity	Mouse
Range	6.25-400pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9R045
Sensitivity	2.45pg/mL
Entrez	26360
Expiration Date	Please enquire.