

Product Datasheet

Rat ACP5 ELISA Kit (orb1200248)

Description	Rat Acid Phosphatase 5, Tartrate Resistant (ACP5) ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P29288
Sensitivity	0.29ng/mL
Entrez	25732
Expiration Date	Please enquire.