

Product Datasheet

Mouse ACE2 ELISA Kit (orb1200252)

Description	Mouse Angiotensin I Converting Enzyme 2 (ACE2) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8R0I0
Sensitivity	24.4pg/mL
Entrez	70008
Expiration Date	Please enquire.