

Product Datasheet

Goat TNFa ELISA Kit (orb1200267)

Description	Goat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Goat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13296
Sensitivity	5.9pg/mL
Entrez	100861232
Expiration Date	Please enquire.