

Product Datasheet

Equine IL1a ELISA Kit (orb1200279)

Description	Equine Interleukin 1 Alpha (IL1a) ELISA Kit
Reactivity	Equine
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q28385
Sensitivity	27pg/mL
Entrez	100064969
Expiration Date	Please enquire.