

Product Datasheet

Equine IL17 ELISA Kit (orb1200282)

Description	Equine Interleukin 17 (IL17) ELISA Kit
Reactivity	Equine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F6U832
Sensitivity	6.0pg/mL
Entrez	100059131
Expiration Date	Please enquire.