

Product Datasheet

Equine IL13 ELISA Kit (orb1200283)

Description	Equine Interleukin 13 (IL13) ELISA Kit
Reactivity	Equine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	B6C802
Sensitivity	6.6pg/mL
Entrez	100034113
Expiration Date	Please enquire.