

Product Datasheet

Goat BMP4 ELISA Kit (orb1200288)

Description	Goat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Goat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	A6N998
Sensitivity	0.148ng/mL
Entrez	100860789
Expiration Date	Please enquire.