

Product Datasheet

Human TGM2 ELISA Kit (orb1200336)

Description	Human Transglutaminase 2, Tissue (TGM2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P21980
Sensitivity	0.055ng/mL
Entrez	7052
Expiration Date	Please enquire.