

Product Datasheet

Human SCCA1 ELISA Kit (orb1200374)

Description	Human Squamous Cell Carcinoma Antigen 1 (SCCA1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29508
Sensitivity	34pg/mL
Entrez	6317
Expiration Date	Please enquire.