

Product Datasheet

Human PTPRJ ELISA Kit (orb1200392)

Description	Human Protein Tyrosine Phosphatase Receptor Type J (PTPRJ) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q12913
Sensitivity	0.117ng/mL
Entrez	5795
Expiration Date	Please enquire.