

Product Datasheet

Human Ntn4 ELISA Kit (orb1200437)

Description	Human Netrin 4 (Ntn4) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9HB63
Sensitivity	32pg/mL
Entrez	59277
Expiration Date	Please enquire.