

Product Datasheet

Human MCSP ELISA Kit (orb1200487)

Description	Human Melanoma Associated Chondroitin Sulfate Proteoglycan (MCSP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q6UVK1
Sensitivity	0.057ng/mL
Entrez	1464
Expiration Date	Please enquire.