

Product Datasheet

Human ITGaM ELISA Kit (orb1200537)

Description	Human Integrin Alpha M (ITGaM) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11215
Sensitivity	0.125ng/mL
Entrez	3684
Expiration Date	Please enquire.