

Product Datasheet

Human ITGa2 ELISA Kit (orb1200540)

Description	Human Integrin Alpha 2 (ITGa2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17301
Sensitivity	0.055ng/mL
Entrez	3673
Expiration Date	Please enquire.